



SN Blood Grouping Teaching Kit

Solution Based

Innovation Speaks

SN Blood Grouping Teaching Kit

Product Code: SNLS -TK217

Number of experiments that can be performed: 100

Duration of Experiment:

Protocol: 15 minutes

Storage Instructions:

- The kit is stable for 12 months from the date of manufacture
- Store Anti A, Anti B and Anti Rh D at 2-8°C
- Other kit contents can be stored at room temperature (15-25°C)

INDEX

Sr. No.	Contents	Page No.
1	Aim	3
2	Introduction	3
3	Principle	3
4	Kit Contents	4
5	Materials Required But Not Provided	4
6	Storage	4
7	Important Instructions	4
8	Procedure	4
9	Observation and Result	5
10	Interpretation	6

Aim:

To determine the blood group and Rh factor of an individual.

Introduction:

Blood grouping is the classification of blood based on the presence or absence of two inherited antigenic substances found on the surface of red blood cells (RBCs). The ABO and Rh systems are the major and most clinically significant blood group systems, and they are considered the most important of all the blood group systems. The ABO blood group system was first discovered by Karl Landsteiner in 1900.

The human ABO blood group system is categorized into four major groups based on the antigens present on the surface of RBCs:

1. "A" group
2. "B" group
3. "AB" group
4. "O" group

Table 1: ABO Blood Group System

Antigens on the surface of Red Blood Cells	Antibodies in the Serum	ABO Blood Group	Genotype
A	Anti B	A	AA or AO
B	Anti A	B	BB or BO
A and B	Neither Anti A nor Anti B	AB	AB
Neither A nor B	Anti A, Anti B, Anti AB	O	OO

The associated Anti A and Anti B antibodies usually belong to IgM class of immunoglobulins.

The Rhesus system (Rh) is the second most important blood group system in humans. The most significant and immunogenic Rhesus antigen is the RhD antigen. The individuals carrying the Rh antigen are considered to have positive blood group whereas those individuals that lack this antigen are considered to have negative blood group.

Principle:

The ABO and Rh blood grouping systems are based on agglutination reactions. When red blood cells (RBCs) carrying specific antigens are exposed to their corresponding antibodies, they interact to form visible agglutination or clumping. In the ABO system, the antigens are O-linked glycoproteins, with the terminal sugar residues on the surface of RBCs determining whether the antigen is A or B. Individuals with blood group A have A antigens on their RBCs and anti-B antibodies in their serum. Similarly, individuals with blood group B have B antigens on their RBCs and anti-A antibodies in their serum. Blood group AB individuals express both A and B antigens on their RBCs and lack anti-A or anti-B antibodies in their serum. In contrast, blood group O individuals do not have A or B antigens on their RBCs but possess both anti-A and anti-B antibodies in their serum. The Rh blood group system involves transmembrane proteins, where the exposed loops on the surface of RBCs interact with corresponding antibodies.

Kit Contents:

SN Blood Grouping Teaching Kit enables rapid identification of ABO blood group and Rh factor depending upon the antigen present on the surface of red blood cells.

Table 2: Enlists the materials provided in this kit with their quantity and recommended storage

Sr. No.	Materials Provided	Quantity	Storage
		100 expts	
1	Anti A Sera	5 ml	2-8°C
2	Anti B Sera	5 ml	2-8°C
3	Anti RhD Sera	5 ml	2-8°C
4	Cavity slide	10 Nos.	RT
5	Disposable Mixing Stick	300 Nos.	RT
6	Blood Lancet	100 Nos.	RT

Materials Required But Not Provided:

Reagents: 70% Alcohol/ Spirit

Other requirements: Cotton

Storage:

SN Blood Grouping Teaching Kit is stable for 12 months from the date of manufacture without showing any reduction in performance. Store the Anti A Sera, Anti B Sera and Anti RhD Sera at 2-8°C. Other contents can be stored at room temperature (15-25°C).


Important Instructions:

1. Before starting the experiment the entire procedure has to be read carefully.
2. Always wear gloves while performing the experiment.
3. Ensure the slide is clean and dry prior to use.
4. Do not allow the antisera reagent dropper to touch the blood sample.
5. The result of the reaction should be interpreted immediately after mixing.
6. Avoid intermixing of the antisera reagents while performing the experiment as it may give false result.

Procedure:

1. Dangle the hand down to increase the flow of blood in the fingers.
2. Clean the fingertip to be pierced with spirit or 70% alcohol (usually ring or middle finger).
3. With the help of the sterile lancet, pierce the fingertip and place one drop of blood in each of the cavities.

4. Add one drop of antiserum into each cavity as shown below:



5. Mix each blood drop and the antiserum using a fresh mixing stick.
6. Observe agglutination in the form of fine red granules within 30 seconds. Anti RhD takes slightly longer time to agglutinate compared to Anti A and Anti B.

Note: Proper care should be taken while disposing the lancet and mixing sticks.

Observation and Result:




Table 3: Determination of blood group and Rh factor based on agglutination seen

Sr. No.	Anti A	Anti B	Anti RhD	Blood Group
Slide 1	✓	×	✓	A +ve
Slide 2	×	✓	✓	B +ve
Slide 3	✓	✓	✓	AB +ve
Slide 4	×	×	✓	O +ve

✓ : Agglutination


× : No agglutination

Interpretation:

- If agglutination is observed when blood is mixed with Anti A reagent, then the individual is said to have blood group "A".
- If agglutination is observed when blood is mixed with Anti B reagent, then the individual is said to have blood group "B".
- If agglutination is observed when blood is mixed with Anti A and Anti B reagent, then the individual is said to have blood group "AB".
- If no agglutination is observed when blood is mixed with Anti A and Anti B reagent, then the individual is said to have blood group "O".
- If agglutination is observed when blood is mixed with Anti RhD reagent, then the individual is said to have "+ve" Rh factor.
- If no agglutination is observed when blood is mixed with Anti RhD reagent, then the individual is said to have "-ve" Rh factor.

Troubleshooting Guide:

Sr.No	Problem	Possible Cause	Solution
1.	False positive result	The antisera reagents mix with each other	Ensure that the antisera reagents are added properly onto the respective cavity without spilling to the sides
		Incubated for a longer time	The results should be read within the time period mentioned in the brochure
2.	No agglutination observed	Antisera not stored under proper conditions	Ensure that the antisera is stored in refrigerator (2-8°C)



Storage temperature



Do not use if package is damaged



**SHANKARANARAYANA
LIFE SCIENCES**
A NEW GENERATION COMPANY



SN Research Unit



SN Manufacturing Unit



SN Diagnostic Unit

Contact Us -

☎ +91- 63613 81129 / 99005 05076

🌐 www.snlifesciences.in

✉ sales@snlifesciences.com / drmondal@snlifesciences.com

🏠 No.57/90, 3rd Road, 4th Phase Bommasandra Industrial Area,
Bengaluru - 560099

2025 - 2026

#startupindia

