




Contact Us -

☎ +91- 80959 61994 / 99005 05076

✉ sales@snlifesciences.com / drmondal@snlifesciences.com

📍 No.57/90, 3rd Road, 4th Phase Bommasandra Industrial Area,
Bengaluru - 560099



#startupindia



SN Plant/Feed RNA extraction Kit

Magnetic Bead Based

SN PLANT/FEED RNA EXTRACTION KIT (Magnetic Bead - Based)

DESCRIPTION:

SN PLANT/FEED RNA Extraction Kit is designed for preparation of high-quality of total RNA from a wide range of PLANT/FEED samples. RNA in the whole homogeneity is selectively absorbed on the Magnetic Bead and other impurities are washed away. This kit provides a simple and efficient way to purify RNA.

CONTENTS OF KIT:

Sl. No	Components	Volume		
		5 Rxn	50 Rxn	100 Rxn
1	RS Buffer	2.5ml	25ml	50ml
2	LS Buffer	3ml	30ml	60ml
3	STB Buffer	2.5ml	25ml	50ml
4	W1 Buffer	1.5ml	15ml	30ml
5	Elution Buffer	250µl	2.5ml	5ml
6	SN Magnetic-Bead	250µl	2.5ml	5ml

NOTE: Preparation for first use after receiving the kit (Add 100% Ethanol to Wash Buffer) (5Rxn: 1ml, 50Rxn:10 ml, and 100Rxn: 20ml) mix well and store the buffer at room temperature.

- ❖ If Proteinase K shipped in lyophilized form, upon receiving resuspend with Nuclease free water (5Rxn:100µl,50Rxn: 1ml, 100Rxn: 2ml)
- ❖ If DNase shipped in lyophilized form, upon receiving resuspend with Nuclease free water (5Rxn:50µl, 50Rxn:500µl, 100Rxn: 1ml)
- ❖ **Proteinase -K and DNase , store at -20°C either liquid form or lyophilized form.**

REQUIRED MATERIALS NOT PROVIDED

- ✓ 100% Ethanol & 70 % Ethanol
- ✓ Dry bath
- ✓ 1.5ml Centrifuge tubes
- ✓ micro centrifuge
- ✓ Vortex

STORAGE / SHIPPING:

- ❖ Shipped at: Ambient Temperature.
- ❖ Storage: All Buffers can be stored at Room temperature.

SPECIFICATIONS:

Principle	Mag-Bead
Recommended Input Amount :200mg of PLANT/FEED sample	
Binding Capacity : ~20-30µg genomic RNA	
Elution Volume : 50µl	
Purity: A260/280 - 1.8±0.1, A260/230 - 2.0±0.1	
Compatible Downstream Applications : PCR, Cloning, Next generation sequencing etc.	
Expected Yield : depending upon sample	

PROCEDURE:

- 1) Weigh 200mg of clean PLANT/FEED sample and transfer to mortar and pestle grind like a fine paste by adding 500µl of RS Buffer. Transfer the above mixture to 1.5 ml centrifuge tube using a clean spatula.
- 2) Add 600µl of LS buffer into the sample along with 20µl of Proteinase K and 10µl of DNase . Vortex for 1min.
- 3) Incubate at 60°C for 15mins with intermittent vortexing for 3-4 times to lysed the cells. Centrifuge the tube for 2mins at 10,000rpm, then transfer the supernatant to a fresh centrifuge tube.
- 4) Add 500µl of 100% ethanol to the supernatant, keep the tube for 2 mins at room temperature, and observe precipitation.
- 5) Transfer the complete precipitate to a fresh centrifuge tube, add 50µl of SN Magnetic bead and 500µl of STB buffer, mix well by vortexing for 30secs then keep at room temperature for 10mins.
- 6) Place the tube upon the magnetic rack (let the bead attracts towards magnet). Then remove the supernatant without disturbing the bead.
- 7) Add 500µl of W1 buffer mix well by vortexing for 30secs then keep at room temperature for 10mins.


Note : (W1 buffer concentration per reaction: W1 buffer- 300µl: 100% ETOH- 200µl)

- 8) **Repeat the step (6)**
- 9) Add 500µl of 70% ethanol, mix well by vortexing for 30secs then keep at room temperature for 5mins.

Repeat the step (6)

- 10) Incubate the tube at 56°C for 10mins (with cap open).
NOTE: Elution should be pre-heated for 10mins before adding into the tube.
- 11) Add 50µl of pre-heated Elution to the above tube, then vortex for 30secs and incubate at 56°C for 4mins (with cap close)
- 12) Vortex for 10secs and then Centrifuge the tube at 15,000rpm for 2mins to Elute pure RNA (store at -20°C)

FLOW CHART:



TROUBLE SHOOTING:

Problems	Possible reasons	Solutions
Low or none recovery of RNA fragment	Weigh 200mg of PLANT/FEED sample	If product is more, then separate it into multiple tubes.
	Elution of RNA fragment is not efficient	Make sure the pH of Elution Buffer or ddH ₂ O is between 7.5-8.0
		Make sure that the elution solution has been observed by the magnetic beads.
Poor Performance in the downstream applications		Preheat the elution solution at 56°C before use.
	Salt residue remains in eluted RNA	
	Ethanol residue remains in eluted RNA	washing step is done with w1 buffer

IMPORTANT NOTES:

- a) Buffer provided in this kit contains irritants. Wear gloves and lab coat when handling these buffers.
- b) Add required volume of ethanol (96-100%) to wash Buffer before use.
- c) Centrifugation steps are done by a micro centrifuge capable of the speed at 11,000 ~ 15,000 rpm.
- d) Fresh TBE electrophoresis buffer is recommended on electrophoresis for experiments requiring high purity.
- e) Then excising the agarose gel, make sure to shorten the time of exposure to ultraviolet irradiation.